

ULTRA FELTAC® and VELCRO® Heavy-Duty Round Felt Pads

Diameter

19 mm 25 mm 48 mm	
-------------------	--

Packaging format

Blister of 12 units	Blister of 24 units	Blister of 48 units
------------------------	------------------------	---------------------

DESCRIPTION

Peel. Stick. Hook.

No mess, no hassle! Thanks to authentic VELCRO® Fastener Brand pads, felt pads are easily changed: just unhook and hook.

ADVANTAGES AND BENEFITS

Slides easily and protects hard surface floors Clean, quick and safe Made in Canada

AVAILABLE PRODUCTS

Product #	Diameter	Packaging format
23231	19 mm	Blister of 48 units
23232	25 mm	Blister of 24 units
23233	48 mm	Blister of 12 units

TECHNICAL SPECIFICATIONS

Recommended Surfaces	Tile/Ceramic/Linoleum, Hardwood/Laminate	
Mounting Type	Self-adhesive	
Shape	Round	
Material	Polyester	
Color / Finish	White, Beige	
Brand	ULTRA FELTAC®	

IMPORTANT INFORMATION

Once the FELTAC® pad is used, simply unhook from the VELCRO® and replace with a new FELTAC® pad. DO NOT turn the FELTAC® pad around. Care must be taken when installing pads to ensure effectiveness; not all furniture sits squarely on floor. Check periodically for wear; felt should be replaced accordingly.

APPLICATION

1. Thoroughly clean the bottom of the leg to remove dirt build-up; leave to dry 2. Remove the protective paper from the VELCRO® (DO NOT touch the adhesive); glue under the leg (CANNOT be repositioned) 3. Apply the FELTAC® pad to the hook side of the VELCRO®



DISCLAIMER

Responsibility is limited to the cost of the product.

PRODUCT PHOTOS

